



Golden Touch
SURGICO

Laryngoscope & Ent Set Catalogue



info@goldentouchsurgico.com
www.goldentouchsurgico.com

Company Profile

Golden Touch Surgico is a family owned company, dedicated to giving you the personalized service you deserve. We have been supplying the medical industry with the finest in Surgical, Dental instrumentation, Laryngoscope, ENT Devices Hospital Furniture and Stainless Steel Hollowware for all specialties for over 40 years. Our office is located in Sialkot, Pakistan, we offer complete in-house repair service and custom manufacturing in addition to our already broad line of products.

When you buy from GTS, you are guaranteed High Quality products. Not only the best products, but the best service as well. If you ever have a question, we're here to help. We will be more than happy to give you the support you need.

Certifications & Standards:

As written on Surgical Instruments Catalogue To
Our Existing Customers

We would like to thank you for continuing to stay with the best. We appreciate your business. To
Prospective Customer

We look forward to starting a lasting relationship with you.

Guarantee

Golden Touch Surgico are guaranteed to be free of defects in materials and workmanship. Any product that is determined to be defective will be either repaired or replaced, at our discretion. Product not properly cared for or used for purposes not intended will not be covered. Golden Touch Surgico liability under this guarantee is limited to repair or replacement of defective items.

OEM Services

We have been providing OEM services to various bands in the health care industry all over the world for a long time

OEM SERVICES PROVIDED, WITH THE UTMOST DISCRETION



GT-500-01
Mc - Intosh Blade Standard Illumination
Size: 0



GT-500-02
Mc - Intosh Blade Standard Illumination
Size: 1



GT-500-03
Mc - Intosh Blade Standard Illumination
Size: 2



GT-500-04
Mc - Intosh Blade Standard Illumination
Size: 3



GT-500-05
Mc - Intosh Blade Standard Illumination
Size: 4

GT-500-06
Handle Standard Illumination
Size: Adult
Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 0, 1, 2 and Adult Handle	▪ GT-500-07
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ GT-500-08
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ GT-500-09
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ GT-500-10

Sets available in hard case or soft pouch

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-11
Mc - Intosh Blade Fiber Optic Illumination
Size: 0



GT-500-12
Mc - Intosh Blade Fiber Optic Illumination
Size: 1



GT-500-13
Mc - Intosh Blade Fiber Optic Illumination
Size: 2



GT-500-14
Mc - Intosh Blade Fiber Optic Illumination
Size: 3



GT-500-15
Mc - Intosh Blade Fiber Optic Illumination
Size: 4

GT-500-16
Handle Fiber Optic Illumination
Size: Adult
Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	GT-500-17
Set of 3 Blades Size 1, 2, 3 and Adult Handle	GT-500-18
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	GT-500-19
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	GT-500-20

Sets available in hard case or soft pouch

Specification:


- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-21
Miller Blade Standard Illumination
Size: 00



GT-500-22
Miller Blade Standard Illumination
Size: 0



GT-500-23
Miller Blade Standard Illumination
Size: 1



GT-500-24
Miller Blade Standard Illumination
Size: 2



GT-500-25
Miller Blade Standard Illumination
Size: 3

GT-500-26
Handle Standard Illumination
Size: Child / Infant
Handle options: Adult & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Child Handle	▪ GT-500-27
▪ Set of 3 Blades Size 0, 1, 2 and Child Handle	▪ GT-500-28
▪ Set of 4 Blades Size 00, 0, 1, 2 and Child Handle	▪ GT-500-29
▪ Set of 4 Blades Size 0, 1, 2, 3 and Child Handle	▪ GT-500-30

Sets available in hard case or soft pouch

Specification:

- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-31
Miller Blade Fiber Optic Illumination
Size: 0



GT-500-32
Miller Blade Fiber Optic Illumination
Size: 1



GT-500-33
Miller Blade Fiber Optic Illumination
Size: 2



GT-500-34
Miller Blade Fiber Optic Illumination
Size: 3



GT-500-35
Miller Blade Fiber Optic Illumination
Size: 4

GT-500-36
Handle Fiber Optic Illumination
Size: Adult
Handle options: Child & Stubby

Standard Sets

Contents	Code
Set of 3 Blades Size 0, 1, 2 and Adult Handle	GT-500-37
Set of 3 Blades Size 1, 2, 3 and Adult Handle	GT-500-38
Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	GT-500-39
Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	GT-500-40

Sets available in hard case or soft pouch

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-47
Handle Fiber Optic Illumination
Size: Adult
Handle options: Child & Stubby

Standard Sets

Contents	Code
▪ Set of 3 Blades Size 00, 0, 1 and Adult Handle	▪ GT-500-48
▪ Set of 3 Blades Size 1, 2, 3 and Adult Handle	▪ GT-500-49
▪ Set of 4 Blades Size 0, 1, 2, 3 and Adult Handle	▪ GT-500-50
▪ Set of 4 Blades Size 1, 2, 3, 4 and Adult Handle	▪ GT-500-51

Sets available in hard case or soft pouch



GT-500-41
Mc - Intosh Blade Fiber Optic Illumination
Size: 00



GT-500-42
Mc - Intosh Blade Fiber Optic Illumination
Size: 0



GT-500-43
Mc - Intosh Blade Fiber Optic Illumination
Size: 1



GT-500-44
Mc - Intosh Blade Fiber Optic Illumination
Size: 2



GT-500-45
Mc - Intosh Blade Fiber Optic Illumination
Size: 3



GT-500-46
Mc - Intosh Blade Fiber Optic Illumination
Size: 4

Specification:

- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-55
Handle Fiber Optic Illumination
Size: Adult



GT-500-52
McCoy Blade Fiber Optic Illumination
Size: 1



GT-500-53
McCoy Blade Fiber Optic Illumination
Size: 2



GT-500-54
McCoy Blade Fiber Optic Illumination
Size: 3

Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	GT-500-56
Set of 2 Blades Size 3, 4 and Adult Handle	GT-500-57
Set of 3 Blades Size 2, 3, 4 and Adult Handle	GT-500-58

Sets available in hard case or soft pouch



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-62
Handle Fiber Optic Illumination
Size: Adult



GT-500-59
McCoy Blade Fiber Optic Illumination
Size: 1



GT-500-60
McCoy Blade Fiber Optic Illumination
Size: 2



GT-500-61
McCoy Blade Fiber Optic Illumination
Size: 3

Standard Sets

Contents	Code
Set of 2 Blades Size 2, 3 and Adult Handle	GT-500-63
Set of 2 Blades Size 3, 4 and Adult Handle	GT-500-64
Set of 3 Blades Size 2, 3, 4 and Adult Handle	GT-500-65

Sets available in hard case or soft pouch



Specification:

- Flexible Tip for difficult Intubation
- Fiber Optic Handle with 3.5mm dia for improved light transmission and longer life.
- Manufactured From AISI 304 Stainless Steel
- All Blades and Handles Compatible with ISO 7376 Standard
- Autoclavable
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available



Head



Handle



Blade



GT-500-66
Brass Otoscope Fiber Optic
Brass Otoscope Conventional



GT-500-67
ENT Set Complete



GT-500-68
Brass Ophthalmoscope

Specification:

- Manufactured From Brass Chrome plated 3x
- Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x C size batteries
- Hard Case and Soft Pouch available
- Fiber Optic Light Source



GT-500-69
Mini Otoscope Fiber Optic
Complete Set of Disposable Cannulas



Color Options

Specification:

- Manufactured From Plastic Material 3x
- Magnifying glass
- Available in Halogen and LED Illumination
- Bulbs 2.5V
- Handles requires 2 x AA size batteries
- Soft Pouch available in matching colors
- Fiber Optic Light Source

Nasal Speculum



GT-500-70



GT-500-71



GT-500-72

Bulbs



GT-500-73



GT-500-74



GT-500-75



GT-500-76

Ear Speculas



GT-500-78



GT-500-79



GT-500-80



GT-500-81

Ear Speculas & Testing Tube



GT-500-82



GT-500-83



GT-500-84



GT-500-85



GT-500-86



GT-500-87



GT-500-88

Colour Option



GT-500-89
Mini Otoscope standard optic



GT-500-90



GT-500-91



GT-500-92

GT

Golden Touch
SURGICO

GOLDEN TOUCH SURGICO.

Sialkot, Pakistan.

+92 300 6148417

Email: info@goldentouchsurgico.com

www.goldentouchsurgico.com



ایوان صنعت و تجارت سialکوت

